



ABSTRACT

The largest current disease-induced loss of vertebrate biodiversity is due to chytridiomycosis and despite the increasing understanding of the pathogenesis, knowledge unravelling the early host-pathogen interactions remains limited. *Batrachochytrium dendrobatidis* (*Bd*)zoospores attach to and invade the amphibian epidermis, with subsequent invasive growth in the host skin. Availability of an *in vitro* assay would facilitate in depth study of this interaction while reducing the number of experimental animals needed. We describe a fluorescent cell-based *in vitro* infection model that reproduces host-*Bd* interactions. Using primary keratinocytes from *Litoria caerulea* and the epithelial cell line A6 from *Xenopus laevis*, we reproduced different stages of host cell infection and intracellular growth of *Bd*, resulting in host cell death, a key event in chytridiomycosis. The presented *in vitro* models may facilitate future mechanistic studies of host susceptibility and pathogen virulence.

EXTERNAL LINK

https://doi.org/10.1371/journal.pone.0225224

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Elin Verbrugghe

Verbrugghe E, Rooij PV, Favoreel H, Martel A, Pasmans F (2019) *In vitro* modeling of *Batrachochytrium dendrobatidis* infection of the amphibian skin. PLoS ONE 14(11): e0225224. doi: 10.1371/journal.pone.0225224

MATERIALS		
NAME ~	CATALOG # V	VENDOR
Distilled Water		
HBSS with calcium and magnesium		
HBSS without calcium and magnesium		
Fetal bovine serum		
Paraformaldehyde	P6148	Sigma Aldrich
Leibovitz's L-15 Medium	11415049	Thermo Fisher
CellTracker™ Red CMTPX Dye	C34552	Thermo Fisher
CellTracker™ Green CMFDA Dye	C7025	Thermo Fisher
Hoechst 33342, Trihydrochloride, Trihydrate, 100 mg	H1399	Thermo Fisher
ProLong™ Glass Antifade Mountant	P36980	Thermo Fisher
Rat Tail Collagen Coating Solution	122-20	Sigma Aldrich

MATERIALS TEXT

24-well tissue culture plates glass coverslips

1 Prepare Cell Medium A:

L15 medium: 70% Distilled water: 20% Fetal bovine serum 10%

2 Prepare Cell Medium B:

L15 medium: 40% Distilled water: 55% Fetal bovine serum: 5%

3 Coat coverslips with Rat tail collagen:

Add glass coverslips in a 24-well tissue culture plate. Coat the glass coverslips at 37° C for 2 hours. Therefore, carefully aspirate the Collagen Coating Solution and add $200 \, \mu$ l per well (so per coverslip). After 2 hours, rinse the coated surface twice with HBSS-

. Coated tissue culture ware may be used immediately or air-dried and stored at 4 °C for up to one week.

4 Staining of Primary Keratinocytes (PAK): (Work in the dark)

- Add 1 ml of 3 µM Celltracker green CMFDA in cell medium A to the PAK
- Resuspend carefully
- Incubate for 30 min at 20°C 5% CO2
- Centrifuge for 5 min at 1500 rpm and 20°C
- Resuspend the pellet in 70% HBSS- and in meantime count the cells
- Centrifuge for 5 min at 1500 rpm
- Resuspend the cells in cell medium A to reach a concentration of 10e5 cells per mL
- Seed the cells at a concentration of 10e5 cells per well (so add 1 mL) which contains a collagen-coated glass coverslip
- Let the PAK attach for 1 hour at 20°C and 5% CO2
- Afther 1 hour, wash the cells with 70% HBSS+
- The cells are now ready to be exposed to Bd spores

IMORTANT NOTE: Include a control well that can be used to check the cells throughout the entire protocol via light microscopy.

5 1

Staining of Bd zoospores: (Work in the dark)

- Isolate Bd zoospores and spin them down for 5 min at 3000 rpm (20°C)
- Remove the supernatant and resuspend them in 3 μ M celltracker red CMPTX solution for 45 min at 20°C. Important here is that the 3 μ M celltracker solution should be diluted in CELL MEDIUM B in order to ensure the motility of the zoospores.
- After 45 min, spin the spores down for 3 min at 3000 rpm
- Resuspend them in cell medium B and count the spores
- Dilute the spores to a concentration of 10e6 spores/ml in cell medium B

IMPORTANT NOTE: Bd zoospores lose their motility when exposed to cell medium A.

IMORTANT NOTE: Include a control well that can be used to check the motility of the spores throughout the entire protocol via light microscopy.

6 Bd infection of PAK cells: (Work in the dark)

- Add 1 ml of the spore suspension (= 10e6 spores) to the wells containing celltracker-labelled PAK cells (on a coverslip). As such the spores are seeded at a MOI of 10:1.
- Incubate for 2 hours at 20°C 5%CO2.
- After 2 hours gently wash the infected cells three times with 70% HBSS+ to remove non-adherent spores
- Replace cell medium B with cell medium A for another 2 hours to assess ADHESION or 22 hours to assess invasion

IMPORTANT: Use the control well to check the adhesion/invasion of the spores to/in the cells via light microscopy.

7 Visualisation of early Bd-PAK interactions: (Work in the dark)

ADHESION: 4 hours post infection (2 hours with cell medium B and 2 hours with cell medium A) INVASION: 24 hours post infection (2 hours with cell medium B and 22 hours with cell medium A)

- Wash the infected cells 3 times with 70% HBSS+
- Fix the infected cells with 0.5 mL of 3% paraformaldehyde for 10 min
- Stain the nulear content with Hoechst
- Wash 2 times with 70% HBSS+
- Mount the coverslips using ProLong Gold antifade mountant
- Use fluorescence microscopy and confocal laser scanning microscopy to analyse Bd-PAK interactions

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited