# Viruses from Hypersaline Waters

### Kenneth M. Stedman, Kate Porter, and Mike L. Dyall-Smith

#### **Abstract**

This is a protocol from:

Stedman, K. M., K. Porter, and M. L. Dyall-Smith. 2010. Chapter 6: The isolation of viruses infecting Archaea. Manual of Aquatic Viral Ecology. Waco, TX:American Society of Limnology and Oceanography. doi:10.4319/mave.2010.978-0-9845591-0-7

Please see the <u>published manuscript</u> for additional information.

Citation: Kenneth M. Stedman, Kate Porter, and Mike L. Dyall-Smith Viruses from Hypersaline Waters. protocols.io

dx.doi.org/10.17504/protocols.io.dph5j5

Published: 26 Feb 2016

#### **Guidelines**

**Authors:** Kenneth M. Stedman<sup>1</sup>, Kate Porter<sup>2</sup>, and Mike L. Dyall-Smith<sup>3</sup>

<sup>1</sup>Department of Biology, Center for Life in Extreme Environments, Portland State University, P.O. Box 751, Portland, OR 97207, USA

<sup>2</sup>Biota Holdings Limited, 10/585 Blackburn Road, Notting Hill Victoria 3168, Australia

<sup>3</sup>Max Planck Institute of Biochemistry, Department of Membrane Biochemistry, Am Klopferspitz 18, 82152 Martinsried, Germany

#### Artificial salt water and medium MGM:

Artificial salt water solutions are designed to mimic the natural concentrated brines where haloarchaea are found. The formulation used by M. Dyall-Smith (described in the online handbook, the *Halohandbook*, <a href="http://www.haloarchaea.com/resources/halohandbook/">http://www.haloarchaea.com/resources/halohandbook/</a>) is based on that described by Rodriguez-Valera et al. (Rodriguez-Valera et al. 1980; Torreblanca et al. 1986). Per liter, it contains:

4 M NaCl

150 mM MgCl<sub>2</sub>

150 mM MgSO<sub>4</sub>

90 mM KCI

3.5 mM CaCl<sub>2</sub>

adjusted to pH 7.5 using ca. 2 mL 1 M Tris-HCl (pH 7.5). At 30% w/v, the total salts are present in a much higher concentration than in seawater, but in approximately the same proportions. Adjustments of Mg<sup>2+</sup>, pH, or other conditions may be necessary for specific haloarchaeal groups. Modified growth medium (MGM) contains 5 g peptone and 1 g yeast extract per liter of salt water. The salt concentration is varied according to the host strain, and is detailed in the <u>Halohandbook</u>.

## **Protocol**