

Introduction to short read assembly

Frank Aylward

Abstract

Here is a short intro on short read assembly using the SPAdes assembler.

Some notes on tool installation:

To install SPAdes or other tools into your PATH you may wish to use a package manager called Miniconda. For this go to the website here and download the 64-bit Linux distribution:

<https://conda.io/miniconda.html>

then, once you have located this file on your computer, type

```
bash Miniconda2-latest-Linux-x86_64.sh
```

and follow the installation instructions. When asked to append the PATH info say "yes".

After this you will need to close your terminal and re-open it before beginning again. Then you should be able to install tools using conda. For the module today we will need to install two tools: the sra-toolkit and an assembler called spades.

```
conda install -c bioconda sra-tools
```

```
conda install -c bioconda spades
```

Citation: Frank Aylward Introduction to short read assembly. **protocols.io**

[dx.doi.org/10.17504/protocols.io.pifdkbn](https://doi.org/10.17504/protocols.io.pifdkbn)

Published: 17 Apr 2018

Protocol

Get the reads

Step 1.

today we will be working with the raw reads from the Staphylococcus phage 812 genome sequencing project.

To get the raw reads we need to use a program called the sra-toolkit. This is a tool maintained by NCBI to allow users to download data easily from the command line.

The toolkit allows users to specify the unique accession number for a project and then download the associated reads.

Here we will also use a command -X, which specifies how many reads we want to download. Since these datasets can be quite large we want to start with a small number. Here we will use 10,000

we also want to use the --split-3 flag, which for Illumina data makes sure the forward and reverse reads are split into separate files.

```
fastq-dump -X 10000 --split-3 SRR6764339
```

Run spades and get the assembly

Step 2.

Now once we have the raw reads we can begin assembling them using a program called SPAdes.
Now actually running the assembly is fairly easy- we just specify the read files, an output folder name, the number of threads we want to use, and the k-mer length.

```
spades.py -1 SRR6764339_1.fastq -2 SRR6764339_2.fastq -o phage -t 4 -k 21 &> log.txt
```


Visualize the assembly

Step 3.

Once you get the assembly you can visualize the De Bruijn graph using a program called Bandage. It can be installed here:

<https://rrwick.github.io/Bandage/>


Once you get it installed you can upload the .fastg file and have Bandage draw the graph. You should get something that looks like this:



Visualize the assembly

Step 4.

You can also color the De Bruijn graph by depth of



coverage on the left-hand side. This can give you an idea of if high-coverage repetitive elements might be causingn problems. This should look like this:

Visualize the assembly

Step 5.

Now go back to the SPAdes command and play around with the k-mer length. How does k-mer length change the topology of the de Bruijn graph?

What about the number of reads used- what if you downloaded 10X fewer or 10X more reads and used those for assembly? How does that change the topology of the de Bruijn graph?