

KAPA Express Extract Reaction Mix

Michael Crone

Abstract

Citation: Michael Crone KAPA Express Extract Reaction Mix. **protocols.io**

dx.doi.org/10.17504/protocols.io.c2pydm

Published: 14 Jun 2015

Protocol

Step 1.

10 X KAPA Express Extract Buffer, **5 µl**

 **AMOUNT**

5 µl Additional info:

Step 2.

1U/µl KAPA Express Extract Enzyme, **1 µl**

 **AMOUNT**

1 µl Additional info:

Step 3.

Water, **44 µl**

 **AMOUNT**

44 µl Additional info:

Step 4.

Sample