

# Thrombin generation assay (CAT)

Sven Van Poucke,Dana Huskens,Kurt Van der Speeten,Mark Roest,Bart Lauwereins,Ming-Hua Zheng,Seppe Dehaene,Joris Penders,Abraham Marcus,Marcus Lancé

## Abstract

TG in plasma was measured with the calibrated automated thrombogram (CAT) assay as developed by Hemker and co-workers [18-20]. Briefly, 80  $\mu$ l platelet poor plasma (PPP) was mixed with 20  $\mu$ l of a mixture containing tissue factor (Dade-Behring) at a final concentration of 1 pM and phospholipid vesicles (f.c. 4  $\mu$ M 20 mol% phosphatidylserine, 60 mol% phosphatidylcholine and 20 mol% phosphatidyl-ethanolamine, Avanti). To calibrator wells, 20  $\mu$ l of calibrator ( $\alpha$ 2macroglobulin- thrombin complex, [19]) was added instead of TF and PL. After 10 minutes of incubation at 37°C, thrombin generation was initiated by the addition of 20  $\mu$ l of the thrombin specific substrate, Z- Gly-Gly-Arg-7-amino-4-methylcoumarin (f.c. 416  $\mu$ M, Bachem) and CaCl<sub>2</sub> (f.c. 16.7 mM). Fluorescence was measured with a Fluoroscan Ascent reader (Thermo Labsystems) and data were analyzed with dedicated software (Thrombinoscope, Stago) [20]. Thrombin generation was expressed based on endogenous thrombin potential (ETP); lagtime (LT); thrombin peak (TP), time-to-thrombin peak (TTP).

**Citation:** Sven Van Poucke,Dana Huskens,Kurt Van der Speeten,Mark Roest,Bart Lauwereins,Ming-Hua Zheng,Seppe Dehaene,Joris Penders,Abraham Marcus,Marcus Lancé Thrombin generation assay (CAT). **protocols.io**  
dx.doi.org/10.17504/protocols.io.ki7cuhn

**Published:** 07 Nov 2017

## Protocol