



May 10,
2019

Working

U Mass - Lipid metabolism [↗](#)

Jason Kim¹

¹University of Massachusetts

[dx.doi.org/10.17504/protocols.io.x3gfqjw](https://doi.org/10.17504/protocols.io.x3gfqjw)

Mouse Metabolic Phenotyping Centers

Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

Summary:

Lipid metabolism is estimated by measuring systemic clearance of [1-¹⁴C] palmitate following a bolus injection in awake mice. Lipid metabolism is altered in obese mice.

EXTERNAL LINK

<https://mmpc.org/shared/document.aspx?id=145&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
Palmitic Acid [1- ¹⁴ C]	NEC075H250UC	Perkin Elmer		

- 1 Survival surgery is performed to establish a chronic indwelling catheter at 5~6 days prior to experiment for intravenous infusion. (refer to M1023: Surgery-jugular vein cannulation)
- 2 Mice are fasted overnight (~15 hours) or for 5 hours prior to the start of experiment.
- 3 Place a mouse in a rat-size restrainer with its tail tape-tethered at one end.
- 4 Expose and flush the intravenous catheter using saline solution. Then, connect the catheter to the CMA Microdialysis infusion pump.
- 5 Collect plasma sample (10 µl) before the start of experiment (basal-0 min) to measure basal glucose levels.
- 6 Administer a bolus intravenous injection of 20 µCi of [1-¹⁴C] palmitate to start the experiment.
- 7 Rapidly collect plasma samples (10 µl each) at 0.5, 1, 2, 3, 4, and 5 min after injection to measure systemic [1-¹⁴C] palmitate concentrations.
- 8 At the end of experiment, mice are euthanized, and tissues may be collected for further studies.

- 9 For data analysis, plasma [1-¹⁴C] palmitate levels vs. time after palmitate injection are plotted, and area-under-curve may be calculated to estimate systemic clearance of labeled-palmitate and lipid metabolism.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited