

# Immunofluorescence assay for detection of ZIKA virus envelope protein

**Elen Mello de Souza, Audrien Alves Andrade, Pedro Paulo de Abreu Manso, Fernanda Cunha Jácome, Marcos Alexandre Nunes da Silva, Gabriela Cardoso Caldas, Ortrud Monika Barth, Debora Ferreira Barreto-Vieira**

## Abstract

**Citation:** Elen Mello de Souza, Audrien Alves Andrade, Pedro Paulo de Abreu Manso, Fernanda Cunha Jácome, Marcos Alexandre Nunes da Silva, Gabriela Cardoso Caldas, Ortrud Monika Barth, Debora Ferreira Barreto-Vieira Immunofluorescence assay for detection of ZIKA virus envelope protein. **protocols.io**

dx.doi.org/10.17504/protocols.io.jqycmxw

**Published:** 05 Sep 2017

## Protocol

### Step 1.

Seed Vero cells on coverslips into 24- well plate at a density of  $2 \times 10^5$  cells/ well.

### Step 2.

Infect with 50  $\mu$ l of the viral suspension at the multiplicity of infection (MOI) 1 for 1h / 37°C to allow for virus adsorption, under agitation every 15 minutes.

### Step 3.

Wash 3 times the cell monolayer with Phosphate-buffered saline (PBSx1) to remove non-adsorbed viral particles.

### Step 4.

Follow the infection up to 72 hours.

### Step 5.

Each 24 h remove two coverslips out the plate.

### Step 6.

Fixation with 500  $\mu$ l of the 4% paraformaldehyde (PFA) for 20 minutes at 4°C.

### Step 7.

Wash 3 times the coverslips with PBS.

### **Step 8.**

Permeabilize the cells with 200 µl of the 0.1% Triton X-100 diluted in PBS for 5 minutes at room temperature.

### **Step 9.**

Wash 3 times the coverslips with PBS.

### **Step 10.**

Block no-specificity binding with 200 µl of the 3% bovine serum albumin (Sigma-Aldrich) for 30 minutes at room temperature.

### **Step 11.**

Wash 3 times the coverslips with PBS.

### **Step 12.**

Incubate with 50 µl of the mouse IgG2a monoclonal antibody anti-ZIKA virus envelope protein (4G2) diluted in PBS (1: 4) for 60 minutes at 37°C.

### **Step 13.**

Wash 3 times the coverslips with PBS.

### **Step 14.**

Incubate with 50 µl of the secondary antibody Alexa Fluor 488 (1: 500) (Invitrogen) for 60 minutes at 37°C.

### **Step 15.**

Wash 3 times the coverslips with PBS.

### **Step 16.**

Incubate with 50 µl of the 10 µg/ ml 4,6-diamidino-2-phenylindole (DAPI) (Sigma-Aldrich) diluted in PBS for 5 minutes at 37°C, for visualization of the cell nucleus.

### **Step 17.**

Mount the coverslips with 10 µl of the 2.5% 1,4-diazabicyclo(2,2,2)-octane (DABCO) (Sigma-Aldrich) to prevent loss of fluorescence.