

WH7803 PCR Protocol

Li Deng

Abstract

Citation: Li Deng WH7803 PCR Protocol. protocols.io

dx.doi.org/10.17504/protocols.io.dak2cv

Published: 21 Jan 2016

Guidelines

Primers:

rpoC1_103: CAC TGC GGC AAA TAC AAA CG

rpoC1_360: TTT CAG GTC CTT GTG GTC AC

Reaction Mix:

Material	Amt.	Final Conc.
Water	14.25µl	
10x Buffer (w/Mg)	2.5µl	=1x
dNTP mix, 2.5mM each	5.0µl	=500 μ M ea.
Primer rpoC1_103, 10μM	0.5μΙ	$=0.2\mu M$
Primer rpoC1_360, 10mM	0.5μΙ	$=0.2\mu M$
La-Taq 5U/μl	0.25µl	$=0.05U/\mu I$
DNA (WH7803)	1.0-2.0µl	
	Total = 25μ l	

Cycling:

95°C	5 min	
95°C	45 sec	
50°C	45 sec	40 cycles
70°C	60 sec	io cycles
70°C	10 min	
15°C	Hold	

Protocol