

1M Citrate and Magnesium Resuspension Buffer

Lauren Chittick

Abstract

Citation: Lauren Chittick 1M Citrate and Magnesium Resuspension Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.exybfpw

Published: 11 May 2016

Protocol

Step 1.

Combine 1.47g sodium citrate dihydrate and 0.05g $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ into 5 ml molecular biology grade water. Add NaOH as necessary to bring pH to 6.5.