



Oct 06, 2019

Virus titration

Norfitriah Mohamed Sohaimi¹, Mohd Hair Bejo¹, Abdul Rahman Omar¹, Aini Ideris¹, Nurulfiza Mat Isa²

¹Faculty of Veterinary Medicine, Universiti Putra Malaysia, Serdang, Selangor, Malaysia., ²Faculty of Biotechnology and Biomolecular Sciences, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.

1 Works for me dx.doi.org/10.17504/protocols.io.7zhhp36

Norfitriah Mohamed Sohaimi 🕢

MATERIALS

NAME VENDOR VENDOR VENDOR Service Culture plates with lids 3596 Corning

MATERIALS TEXT

Viral supernatant from passage 5 (CEL5), CEL10, CEL15, CEL20, CEL25, CEL30 and CEL35 were determined for titration by 50% tissue culture infective dose ($TCID_{50}$) as described by Reed and Muench [1938]. Supernatant was diluted for 10 folds serial dilution and inoculated onto fresh CEL cells in 96 wells plate and monitored daily for CPE. Virus titre was calculated and expressed as $TCID_{50}/mL$.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited