

Digestion Mixture for M0302

New England Biolabs

Abstract

Citation: New England Biolabs Digestion Mixture for M0302. [protocols.io](https://doi.org/10.17504/protocols.io.cqrvv5)
[dx.doi.org/10.17504/protocols.io.cqrvv5](https://doi.org/10.17504/protocols.io.cqrvv5)
Published: 05 Feb 2015

Protocol

Step 1.

DNA, **200 ng**

Step 2.

10X NEBuffer 2, **2 µl** (1X)

 **AMOUNT**

2 µl Additional info:

Step 3.

Nuclease-free Water to **19 µl**