

Dishwashing

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt Dishwashing. protocols.io

dx.doi.org/10.17504/protocols.io.e7xbhpn

Published: 21 Jun 2016

Protocol

Step 1.

Rinse all items with MilliQ before placing in dishwasher

Step 2.

Use appropriate dishwashing rack for items. Make sure black hub is lined up in dishwasher

Step 3.

Add a scoop of dishwashing detergent from white bucket under sink to detergent tray and close tray

Step 4.

Close dishwasher door securely. Turn dial to appropriate wash cycle and press start.