



Microbiome assay 👄

PLOS Genetics

Wen Aw1

¹z3314717@unsw.edu.au

dx.doi.org/10.17504/protocols.io.shjeb4n



Wen Aw



EXTERNAL LINK

https://doi.org/10.1371/journal.pgen.1007735

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Aw WC, Towarnicki SG, Melvin RG, Youngson NA, Garvin MR, Hu Y, Nielsen S, Thomas T, Pickford R, Bustamante S, Vila-Sanjurjo A, Smyth GK, Ballard JWO (2018) Genotype to phenotype: Diet-by-mitochondrial DNA haplotype interactions drive metabolic flexibility and organismal fitness. PLoS Genet 14(11): e1007735. doi: 10.1371/journal.pgen.1007735

PROTOCOL STATUS

Working

- Extract total DNA from 6 sets of 4 guts of early third instar larvae using the NucleoSpin® Tissue XS DNA extraction kit (Machery-Nagel, Düren, Germany) following the manufacturer's instructions.
- Amplify the V4 region of the 16S rRNA gene by PCR in duplicate for each sample, using the 515 forward 5'-GTGCCAGCMGCCGCGGTAA-3' and 806 reverse 3'-GGACTACHVGGGTWTCTAAT-5' primer.
- Pool duplicate PCR products and sequence using a MiSeq platform with 2 x 250 bp chemistry 3
- Merge paired-end sequences into contigs, and cluster into Operational Taxonomic Units (OTUs) at 97% sequence similiarty using Mothur (http://www.mothur.org) and its standard operating procedure for MiSeq data, but with minor changes (here, singleton contigs were removed after the pre-clustering step).
- Taxonomically classify OTUs using the Ribosomal Database Project taxonomy (http://rdp.cme.msu.edu) with 60% confidence 5 threshold, and compare communities at a genus level.
- Rarefy each sample to the same total number of sequences (n=41,859) to account for differenti sequencing depth, then convert to relative abundance for analysis.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

✓ protocols.io