

# Digestion Solution

Kristin Corrier/Nathan VerBerkmoes

## Abstract

For use in the "[FASP Kit Protocol-ORNL Developed for Bacteriophage](#)"

**Citation:** Kristin Corrier/Nathan VerBerkmoes Digestion Solution. **protocols.io**  
dx.doi.org/10.17504/protocols.io.dek3cv

**Published:** 21 Jan 2016

## Protocol

### Step 1.

75uL ammonium bicarbonate solution added to 20µg trypsin (1 tube).