



Jun 12,
2019

Bold SW 1 NV modified medium

Roscoff Culture Collection¹

¹Station Biologique, CNRS-Sorbonne Université, Roscoff France

Working

dx.doi.org/10.17504/protocols.io.s23eggn

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France



ABSTRACT

Medium used for some eukaryotes

Adapted from CCAP web site : <https://www.ccap.ac.uk/media/documents/BB.pdf>

See also : <http://www.ife.ac.uk/ccap/MediaRecipes.html#bb>

MATERIALS

NAME

Red Sea Salt

CATALOG #

View

VENDOR

Red Sea

1 Change from standard Bold medium

Volume for 500 ml	Component	Concentration Stock Solutions		
5 ml	CaCl ₂	CaCl ₂ ·2H ₂ O	0.5 gr	200 ml
5 ml	NaNO ₃	NaNO ₃	16 gr	200 ml
5 ml	MgSO ₄	MgSO ₄ ·7H ₂ O	10 gr	200 ml
10 ml	PO ₄	From PCRS11 media		
5 ml	NH ₄ Cl (K)	From K media		
3 ml	K trace metal	From K media		
50 µl	Vitamines			
467 ml	SW sterile			



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited