

CviRI Buffer C

DAVID DUNIGAN AND IRINA AGARKOVA

Abstract

For use in [CviRI Purification From XZ-6E Virus Infected NC64A Chlorella](#).

Citation: DAVID DUNIGAN AND IRINA AGARKOVA CviRI Buffer C. **protocols.io**

dx.doi.org/10.17504/protocols.io.esybefw

Published: 13 Jun 2016

Protocol

Step 1.

10 mM KHPO₄, pH 7.4

Step 2.

0.5 mM EDTA

Step 3.

7 mM 2-ME, 10% Glycerol