

SDS-PAGE Cracking Buffer

David Dunigan and Irina Agarkova

Abstract

Recipe to make 10ml of 2x stock

Citation: David Dunigan and Irina Agarkova SDS-PAGE Cracking Buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.erpb5n

Published: 13 Jun 2016

Guidelines

1M Tris-HCl pH 6.8	1.0 ml
SDS	0.4 g
Glycerol	2.0 ml
0.5 M EDTA	0.5 ml
Bromophenol Blue	4 mg
β-mercaptoethanol (14.7M)	0.2 ml

* It is better to add BME right before using buffer.

Protocol