

May 08, 2019

Working

## Case - Measurement of Food Consumption and Body Weight

Henri Brunengraber<sup>1</sup>

<sup>1</sup>Case Western Reserve University

dx.doi.org/10.17504/protocols.io.yeqftdw

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org

Lili Liang

ABSTRACT

## Summary:

During the course of development or during a diet study the amount of food and body weights are measured 1-2 times per week.

**EXTERNAL LINK** 

https://mmpc.org/shared/document.aspx?id=269&docType=Protocol

- Mice are marked by ear notches. They are weighed and the weight is recorded.
- 2 At 9:00AM, mice are placed in a clean cage (2-3 mice/cage because singly housed mice are stressed). The amount of food is weighed and recorded.
- 3 The next day at 9:00AM the food and mice are weighed and measurements are recorded.
- 1 To calculate food intake. Subtract day 3 values from day 2 values and divide by the total weight of the mice in the cage.
- 5 The data is represented as an average over the course of diet study or each week over time.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited