OPEN ACCESS



# Immunoblot based assay for simultaneous densitometric determination of ubiquitin forms in Drosophila melanogaster

Ágota Nagy,Levente Kovács,Zoltán Lipinszki,Margit Pál,Péter Deák

# **Abstract**

This protocol describes an immunoassay, originally developed to quantitate ubiquitin in mouse tissues (Oh et al., Anal. Biochem. 443, 153–155), which we adapted to *Drosophila melanogaster*. The method is suitable for the simultaneous determination of total, free and conjugated ubiquitin forms from whole protein extracts by densitometric analysis of Western blots. In this assay, endogenous deubiquitylating enzymes, DUBs, present in the lysates process all conjugated ubiquitins to monoubiquitins, therefore the total ubiquitin content of cell lysates can be determined in the form of monoubiquitins. The free monoubiquitin fraction in turn is determined from similar lysates, but supplemented with a potent DUB inhibitor, NEM. Appropriate samples of these lysates are immunoblotted together with ubiquitin standards that allow the quantification of the different ubiquitin forms by densitometric analysis of the 8,5 kDa monoubiquitin band.

**Citation:** Ágota Nagy, Levente Kovács, Zoltán Lipinszki, Margit Pál, Péter Deák Immunoblot based assay for simultaneous densitometric determination of ubiquitin forms in Drosophila melanogaster. **protocols.io** 

dx.doi.org/10.17504/protocols.io.s5teg6n

Published: 31 Aug 2018

#### **Guidelines**

Buffers and reagents
Tris base
NaCl
Glycine
Tween-20
EDTA
N-Ethylmaleimide (NEM, Sigma-Aldrich)
MG132 (Calbiochem)
EDTA-Free Complete Protease Inhibitor Cocktail (Roche)
DTT

Glycerol

Bromophenol blue

2-mercaptoethanol

Purified ubiquitin (Sigma-Aldrich)

PVDF membrane (Merck Immobilon-P)

Tris-Glycine SDS-PAGE gel

30-40% Acrylamid solution (29:1 acrylamide/bis-acrylamid)

Non-fat milk powder

anti-Ub antibody (rabbit, Dako)

Bovine serum albumin

HRP conjugated goat anti-rabbit secondary antibody

Immobilon Western Chemiluminescent HRP Substrate (Merck)

X-ray films (Fujifilm)

Tetenal X-ray film processing developer and fixative solutions

# **Before start**

Precool the centrifuge to 4 °C. Buffers F and T should be prepared prior to use and kept on ice.

Buffer F:

100 mM Tris, pH 7.6

150 mM NaCl

1 mM EDTA

10 mM N-Ethylmaleimide (NEM)

20 μM MG132

1× EDTA-Free Complete Protease Inhibitor Cocktail

Buffer T:

100 mM Tris, pH 7.6

150 mM NaCl

20 μM MG132 (Calbiochem)

1× EDTA-Free Complete Protease Inhibitor Cocktail (Roche)

2 mM DTT

4 x Laemmli sample buffer:

40% glycerol

240 mM Tris-Hcl pH=6.8

8% SDS

0.04% Bromophenol blue

5 % 2-mercaptoethanol

TBS:

10 mM Tris pH 8.0

150 mM NaCl

#### **Protocol**

# Sample preparation

# Step 1.

Collect 5-6 mg of synchronized animal or tissue samples from *Drosophila melanogaster* in pre-chilled 1.5 ml microfuge tubes.



**AMOUNT** 

5 mg:

# Sample preparation

# Step 2.

Add 100 µl buffer F (in which conjugated ubiquitins remain intact) and 100 µl buffer T (in which all conjugated ubiquitins are converted to monomers) to the samples and homogenize them by plastic tissue grinders.

AMOUNT  $100 \mu l$ : Sample preparation Step 3. Centrifuge the samples at 4 °C, 13000 RCF. **O** DURATION 00:10:00: Sample preparation Step 4. Collect supernatants to new 1.5 ml microfuge tubes. AMOUNT 70 μl: Sample preparation Step 5. Centrifuge the collected samples again at 4 °C, 13000 RCF. **O DURATION** 00:10:00: Sample preparation Step 6. Collect supernatants into new 1.5 ml microfuge tubes. **■** AMOUNT 61 µl: Sample preparation Step 7. Use 1 µl extract to determine the total protein content of the extracts in buffer F. Sample preparation Step 8. Add 4 x Laemmli sample buffer to the protein extracts in buffer F and boil. **■** AMOUNT

20 μl: **O DURATION** 00:05:00:

# Sample preparation

## Step 9.

Incubate the protein extracts in buffer T at 25  $^{\circ}$ 

**O DURATION** 

03:00:00:

# Sample preparation

# Step 10.

Use 1 µl extract to determine the total protein content of the extracts in buffer T.

## Sample preparation

#### Step 11.

Add 4 x Laemmli sample buffer to the protein extracts in buffer T and boil.

**■** AMOUNT

20 μl:

O DURATION

00:05:00:

# Western blot

#### Step 12.

Load 10  $\mu$ l of appropriately diluted total protein extract onto 14 %, 1 mm thick Tris-Glycine SDS-PAGE gel, together with ubiquitin standards of 0.5, 1, 2 and 3 pmols.

#### NOTES

The appropriate loading concentration (determined in Western blot test runs) is the one in which the band intensity of the monoubiquitin appears to be within the band intensity range of the ubiquitin standards.

#### Western blot

#### Step 13.

Perform the SDS-PAGE in a buffer containing 25 mM Tris pH 8.3, 192 mM Glycine and 0.1% SDS.

# Western blot

# Step 14.

Perform standard wet transfer of the proteins onto a PVDF membrane in a blotting apparatus (300 mA constant current, in transfer buffer containing 20 mM Tris pH 8.0, 150 mM Glycine and 20 % Methanol).

**O** DURATION

02:30:00:

#### Western blot

# Step 15.

Block the membrane by rinsing it in 5% non-fat milk in TBS at room temperature.

**O** DURATION

00:30:00:

# Western blot

# Step 16.

Incubate the membrane in anti-ubiquitin antibody, diluted 1:1000 in TBST-B (TBS supplemented with 0.05 % Tween-20 and 1 % bovine serum albumin) at room temperature.

**O DURATION** 

01:00:00:

#### Western blot

## Step 17.

Wash the membrane three times by vigorous shaking in TBS-T.

© DURATION 00:10:00 : 3x

# Western blot

## Step 18.

Incubate the membrane with the secondary antibody, a horseradish peroxidase (HRP) conjugated goat anti-rabbit antibody diluted 1:30000 in TBST-B at room temperature.

© DURATION

01:00:00:

#### Western blot

#### Step 19.

Wash the membrane three times by vigorous shaking in TBS-T.

© DURATION 00:10:00 : 3x

# Western blot

# Step 20.

Incubate the membrane in Immobilon Western Chemiluminescent HRP Substrate for 3 minutes, then cover the membrane with Saran Wrap $^{\text{TM}}$ .

# Western blot

#### Step 21.

In a dark room position a sheet of X-ray film over the membrane for the appropriate exposure time (20-30 seconds), then develop the film immediately using Tetenal X-ray film processing developer and fixative solutions.

#### Data analysis

# Step 22.

Digitalize the developed X-ray films by using a digital camera or a high resolution gel documentation system (such as BioDoc-It<sup>™</sup> 220 Imaging System).

# Data analysis

# Step 23.

Use the Gel Analysis tool of the ImageJ 1.49v software (NIH, Bethesda, Maryland) to determine the density of the bands corresponding to the monoubiquitin and ubiquitin standards. Copy the band intensity data to an excel table for subsequent statistical analysis.

**ELINK:** 

https://imagej.nih.gov/ij/index.html

#### Data analysis

# Step 24.

Use the band intensity data of the ubiquitin standards (y) and their concentration (x) to generate a calibration curve in MS Excel with the XLSTAT, a statistical software and data analysis add-on for Excel. Create a regression line equation by applying the four parameter curve fit model (which can be found at the dose toolbar), then use it to calculate the ubiquitin concentration of the sample of interest.

@ LINK:

https://www.xlstat.com/en/