# **SOC Media**

## **New England Biolabs**

#### **Abstract**

Citation: New England Biolabs SOC Media. protocols.io

dx.doi.org/10.17504/protocols.io.crsv6d

Published: 03 Feb 2015

## **Protocol**

#### Step 1.

**SOB Media** 



## . SOB Media

**CONTACT: New England Biolabs** 

Step 1.1.

2% tryptone

Step 1.2.

0.5% yeast extract

Step 1.3.

10 mM NaCl

Step 1.4.

2.5 mM KCl

Step 1.5.

10 mM MgCl2

Step 1.6.

10 mM MgSO4

#### Step 2.

20 mM glucose