

# DTT solution

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova DTT solution. **protocols.io**

dx.doi.org/10.17504/protocols.io.erabd2e

**Published:** 13 Jun 2016

## Guidelines

Components:

2.25 M DTT, 40 mM KOAc, pH 6.0.

## Protocol

### Step 1.

Add components together and adjust pH as necessary.

### Step 2.

Filter sterilize.

### Step 3.

Store frozen at -20°C.